

**MW- 47 (A) FULLY MOTORISED BED**


Overall Dimension : 2230mm (L) x 1050mm (W)  
 Suitable for mattress : 1980mm (L) x 880mm (W)

Minimum Height : 450mm (Without mattress)  
 Maximum Height : 7500mm (Without mattress)

**Features:**

- 60mm x 30mm CRCA tubular frame work.
- Four section perforated CRCA sheet top.
- Height Adjustment, Backrest, Kneerest, Trendelenburg & R T, Electrically Operated by Handset.
- CPR Function, NCP Box & Battery Back-up.
- Removable Polymer Head and Foot Board.
- Split type Polymer Head and Foot Boards.
- Tuck away type Polymer Railing (set of 4).
- S. S. Telescopic Saline Rod with two locations.
- Central & Directional locking castors.
- Four side corner buffers.
- Urine bag holder.
- Pretreated and Epoxy powder coated.

**• OPTIONAL ACCESSORIES AT EXTRA COST.**

- Infusion pump rod with telescopic I. V. Rod.
- PU Foam Mattress.
- Bed extension facility.



**MW - 47 (B) SEMI MOTORIZED BED**


Overall Dimension : 2230mm (L) x 1050mm (W)  
 Suitable for mattress : 1980mm (L) x 880mm (W)

Minimum Height : 450mm (Without mattress)  
 Maximum Height : 800mm (Without mattress)

**Features:**

- 60mm x 30mm CRCA tubular frame work.
- Four section perforated CRCA sheet top.
- Backrest, Kneerest & Height adjustment are Electrically Operated with Handset
- Removable Polymer Head and Foot Board.
- Tuck away type Polymer Rails (set of 4)
- S. S. Telescopic Saline Rod with two locations.
- Four heavy duty swivel castors of 125mm dia (two with brake).
- Pretreated and epoxy powder coated.

**• OPTIONAL ACCESSORIES AT EXTRA COST.**

- Nursing Control Panel.
- CPR function.
- Trendelenburg & Reverse Trendelenburg.
- Central & Directional locking castors.
- Infusion pump rod with telescopic I. V. Rod.
- PU Foam Mattress.
- Urine bag holder.
- Bed extension facility.



# ACUTE CARE (DELUXE)

## MW - 08 (E) ACUTE CARE DELUXE HEIGHT ADJUSTABLE ICU BED WITH SPLIT TYPE SIDE RAILINGS



Overall Dimension : 2350mm (L) x 1050mm (W) x Ht. adj. from 550mm to 750mm.  
Suitable for mattress : 1980mm (L) x 890mm (W)

### Features:

- 60mm x 30mm CRCA tubular frame work.
- Four section perforated CRCA sheet top.
- Height Adjustment, Backrest, Kneerest, Trendelenburg & Reverse Trendelenburg operated by screw and lever mechanism with four separate S. S. folding handles.
- Polymer moulded Head and Foot Board.
- Split type side Railings (set of 4).
- Tuckaway type Polymer moulded side rails (set of 4)
- S. S. Telescopic Saline Rod with two locations.
- Four heavy duty swivel castors of 125mm dia (two with brake).
- Four side corner buffers.
- Pretreated and epoxy powder coated.

### • OPTIONAL ACCESSORIES AT EXTRA COST.

- Central & Directional Locking Castors.
- Infusion pump Rod with telescopic I.V. Rod.
- Four section PU foam Mattress.
- Urine Bag Holder.



# ACUTE CARE (DELUXE)

## MW - 08 (F) ACUTE CARE DELUXE HEIGHT ADJUSTABLE ICU BED WITH POLYMER RAILINGS (SET OF TWO)



Overall Size : 2350mm (L) x 1050mm (W) x Ht. adj. from 550mm to 750mm.  
Suitable size for mattress : 1980mm (L) x 890mm (W)

### Features:

- 60mm x 30mm CRCA tubular frame work.
- Four section perforated CRCA sheet top.
- Backrest, Kneerest, Trendelenburg, Reverse Trendelenburg & Height Adjustment operated by screw and lever mechanism with four separate S. S. folding handle.
- Polymer moulded Head and Foot Board.
- A pair of new design side rails (1320mm length)
- S. S. Telescopic Saline Rod with four locations.
- Four heavy duty swivel castors of 125mm dia (two with brake).
- Four side corner buffers.
- Pretreated and epoxy powder coated.

### OPTIONAL ACCESSORIES AT EXTRA COST.

- Central & Directional Locking Castors.
- Infusion pump rod with telescopic I.V. Rod.
- Four section PU foam Mattress.
- Urine bag holder.



# ACUTE CARE (DELUXE)

## MW- 08 (D) ACUTE CARE DELUXE HEIGHT ADJUSTABLE ICU BED WITH NEW DESIGN COLLAPSIBLE SIDE RAILINGS



Overall Size : 2350mm (L) x 1050mm (W) x Ht. adj. from 550mm to 750mm.  
Suitable size for mattress : 1980mm (L) x 890mm (W)

### Features:

- 60mm x 30mm CRCA tubular frame work.
- Four section perforated CRCA sheet top.
- Backrest, Kneerest, Trendelenburg, Reverse Trendelenburg & Height Adjustment operated by screw and lever mechanism with four separate S. S. folding handle.
- Polymer moulded Head and Foot Board.
- A pair of new design side rails (1320mm long)
- S. S. Telescopic Saline Rod with four locations.
- Four heavy duty swivel castors of 125mm dia (two with brake).
- Corner Buffers.
- Pretreated and epoxy powder coated.

### • OPTIONAL ACCESSORIES AT EXTRA COST.

- Central & Directional Locking Castors.
- Infusion pump rod with telescopic I.V. Rod.
- Four section PU foam Mattress.
- Urine bag holder.



# ACUTE CARE (DELUXE)

## MW - 07 (D) ACUTE CARE DELUXE FIXED HEIGHT ICU BED WITH NEW DESIGN COLLAPSIBLE SIDE RAILINGS



Overall Size : 2350mm (L) x 1050mm (W) x 550mm (H).  
Suitable size for mattress : 1980mm (L) x 890mm (W)

### Features:

- 60mm x 30mm CRCA tubular frame work.
  - Four section perforated CRCA sheet top.  
Backrest, Kneerest, Trendelenburg & Reverse Trendelenburg operated by screw and lever mechanism with three separate S. S. folding handles.
  - Polymer Head and Foot Board.
  - A pair of new design side railings (1320mm length).
  - S. S. Telescopic Saline Rod with four locations.
  - Four heavy duty swivel castors of 125mm dia (two with brake).
  - Four side corner buffers.
  - Pretreated and epoxy powder coated.
- **OPTIONAL ACCESSORIES AT EXTRA COST.**
  - Central & Directional Locking Castors.
  - Infusion pump rod with telescopic I.V. Rod.
  - Four section PU foam Mattress.
  - Urine bag holder.



# ACUTE CARE (RECOVERY)

## MW - 05 (A) ACUTE CARE RECOVERY BED WITHOUT RAILINGS & CASTORS



## MW - 05 (B) ACUTE CARE RECOVERY BED WITH SWING TYPE SIDE RAILINGS & CASTORS



Overall Size 1) : 2250mm (L) x 1020mm (W) x 590mm (H).  
Suitable size for mattress : 1905mm (L) x 890mm (W)

Overall Size 2) : 2170mm (L) x 930mm (W) x 590mm (H).  
Suitable size for mattress : 1830mm (L) x 810mm (W)

### Features:

- 60mm x 30mm CRCA tubular frame work.
- Two section perforated CRCA sheet top.
- Backrest, Trendelenburg & Reverse Trendelenburg operated by screw and lever mechanism with attached S. S. folding handles.
- 32mm Dia. Stainless Steel bows with coloured metal panel.
- A pair of swing away side railings.
- Four locations for I.V. Rod.
- Four heavy duty swivel castor of 125mm dia (two with brake).
- Pretreated and epoxy powder coated.

- OPTIONAL ACCESSORIES AT EXTRA COST.
- Saline rod.
- PU foam mattress.



**MW- 08 (B): ACUTE CARE  
HEIGHT ADJUSTABLE ICU BED WITH FULL LENGTH COLLAPSIBLE SIDE RAILINGS**

Overall Size1 :  
2350mm (L) x 1020mm (W) x Ht. Adj. from 550mm to 750mm.  
Suitable size for mattress : 2250mm (L) x 890mm (W)

**Features:**

- 60mm x 30mm CRCA tubular frame work.
- Height adjustment, Backrest, Kneerest Trendelenburg and Reverse Trendelenburg operated by screw and lever mechanism with four separate attached S.S. folding handles.
- 32mm dia stainless steel bows with coloured metal panel
- Four section perforated CRCA sheet top
- A pair of Full Length Collapsible side rails.
- Four heavy duty swivel castors of 125mm dia (two with brake).
- S. S. Telescopic Saline Rod with four locations.
- Four side Corner Buffers.
- Pretreated and epoxy powder coated.

Overall Size2 :  
2250mm (L) x 930mm (W) x Ht. Adj. from 550mm to 750mm.  
Suitable size for mattress : 1905mm (L) x 810mm (W)

**OPTIONAL ACCESSORIES AT EXTRA COST.**

- Infusion pump Rod with telescopic I.V. Rod.
- Four Section PU foam mattress.
- Urine bag holder



**MW- 08 ( C ): ACUTE CARE**  
**HEIGHT ADJUSTABLE ICU BED WITH NEW DESIGN COLLAPSIBLE SIDE RAILINGS**



Overall Size1 :  
 2350mm (L) x 1020mm (W) x Ht. Adj. from 550mm to 750mm.  
 Suitable size for mattress : 1980mm (L) x 890mm (W)

**Features:**

- 60mm x 30mm CRCA tubular frame work.
- Four section perforated CRCA sheet top
- Height adjustment, Backrest, Kneerest Trendelenburg, Reverse Trendelenburg operated by screw and lever mechanism with four separate attached S.S. folding handles.
- 32mm dia stainless steel bows with coloured metal panel.
- New Design collapsible Railings (1320 mm long).
- S.S. Telescopic Saline Rod with four locations.
- Four heavy duty swivel castors of 125mm dia.
- Four side Corner Buffers.
- Pretreated and epoxy powder coated.

Overall Size2 :  
 2250mm (L) x 930mm (W) x Ht. Adj. from 550mm to 750mm  
 Suitable size for mattress : 1905mm (L) x 810mm (W)

**OPTIONAL ACCESSORIES AT EXTRA COST.**

- Infusion pump rod with attachment.
- PU foam Mattress.




**MW - 07 (A): ACUTE CARE BED FIXED HEIGHT WITH SWING AWAY RAILINGS**

Overall Size1 :

2350mm (L) x 1020mm (W) x Ht. Adj. from 590mm (H)

Suitable size for mattress : 1985mm (L) x 890mm (W)

**Features:**

- 60mm x 30mm CRCA tubular frame work.
- Backrest, Kneerest, Trendelenburg, Reverse Trendelenburg operated by screw and lever mechanism with attached S.S. folding handles.
- 32mm dia stainless steel bows with coloured metal panel
- S.S. Telescopic Saline Rod with four locations.
- Four Side Corner Buffers.

Overall Size2 :

2250mm (L) x 930mm (W) x Ht. Adj. from 590mm(H)

Suitable size for mattress : 1905mm (L) x 810mm (W)

- A pair of Full Length Collapsible side railings.
- Four heavy duty swivel castors of 125mm dia.
- Pretreated and epoxy powder coated.

**OPTIONAL ACCESSORIES AT EXTRA COST.**

- Infusion pump rod with telescopic I.V. Rod.
- PU foam Mattress.
- Urine Bag Holder.


**VS - MI - 07 (B): ACUTE CARE BED FIXED HEIGHT WITH FULL LENGTH COLLAPSIBLE RAILINGS**

Overall Size1 :

2350mm (L) x 1020mm (W) x 590mm (H).

Suitable size for mattress : 1980mm (L) x 890mm (W)

**Features:**

- 60mm x 30mm CRCA tubular frame work.
- Backrest, Kneerest, Trendelenburg, Reverse Trendelenburg operated by screw and lever mechanism with three separate S.S. folding handles.
- S.S. Telescopic Saline Rod with four locations.
- A pair of Full Length Collapsible side rails
- Pretreated and epoxy powder coated.

- Four section perforated CRCA sheet top.
- 32mm dia stainless steel bows with coloured metal panel.
- A pair of Full Length Collapsible side rails.
- Four heavy duty swivel castors of 125mm dia. (Two with brake).

**OPTIONAL ACCESSORIES AT EXTRA COST.**

- Infusion pump rod with telescopic I.V. Rod.
- PU foam Mattress.
- Urine bag holder.



**MW - 07 (C) : ACUTE CARE FIXED HEIGHT BED  
WITH NEW DESIGN COLLAPSIBLE SIDE RAILINGS**

Overall Size1 :  
2350mm (L) x 1020mm (W) x 590mm (H).  
Suitable size for mattress : 1980mm (L) x 890mm (W)  
Overall Size2 :  
2250mm (L) x 930mm (W) x 590mm (H).  
Suitable size for mattress : 1905mm (L) x 810mm (W)

**Features:**

- 60mm x 30mm CRCA tubular frame work.
- Four section perforated CRCA sheet top
- Backrest, Kneerest, Trendelenburg & Reverse Trendelenburg operated by screw and lever mechanism with three separate attached S.S. folding handles.
- A pair of New Design Collapsible side railings (1320mm length).
- S.S. Telescopic Saline Rod with four locations.
- Four heavy duty swivel castors of 125mm dia (two with brake).
- Four side Corner Buffers.
- Pretreated and epoxy powder coated.

**OPTIONAL ACCESSORIES AT EXTRA COST.**

- Infusion pump Rod with telescopic I.V. Rod.
- PU foam Mattress.
- Urine bag holder.



**MW - 03 (C): SEMI FOWLER BED (DELUXE)-PRIMA****Overall Size1 :**

2060mm (L) x 1050mm (W) x 590mm (H).

Suitable size for mattress : 1905mm (L) x 890mm (W)

**Features:**

- 60mm x 30mm CRCA tubular frame work.
- Two section perforated CRCA sheet top.
- Backrest operated by screw and lever mechanism with attached S.S. folding handle.
- Polymer moulded head & foot board
- A pair of new design collapsible rails (1320mm length)
- Four heavy duty swivel castors of 125mm Dia. (two with brake).
- Four side corner buffers.
- Four locations for I.V. Rod.
- Pretreated and epoxy powder coated.

**OPTIONAL ACCESSORIES:**

- a) Saline rod.
- b) Urine bag holder.
- c) PU foam mattress.



**MW - 04 (B) FOWLER BED (DELUXE) - PRIMA****Overall Size1 :**

2350mm (L) x 1050mm (W) x 590mm (H).

Suitable size for mattress : 1980mm (L) x 890mm (W)

**Features:**

- 60mm x 30mm CRCA tubular frame work.
- Four section perforated CRCA sheet top.
- Backrest and Kneerest operated by screw and lever mechanism with s.s. folding handle.
- Polymer moulded head and foot boards.
- A pair of new design Colapsible side railings (1320mm length)

- Four heavy duty swivel castors of 125mm Dia, (two with brake)
- Four locations for I.V. Rod.
- Pretreated and epoxy powder coated.

**OPTIONAL ACCESSORIES:**

- a) Saline Rod.
- b) PU foam mattress.





### MW - 04 (A) FOWLER BED (STANDARD)

Overall Size1 :

2160mm (L) x 925mm (W) x 590mm (H).

Suitable size for mattress : 1980mm (L) x 890mm (W)

#### Features:

- 60mm x 30mm CRCA tubular frame work.
- Four section perforated CRCA sheet top.
- Backrest and Kneerest operated by screw and lever mechanism with s.s. folding handle.
- 32mm dia stainless steel bows with coloured metal panel.
- Four locations for I.V. Rod.
- Pretreated and epoxy powder coated.



Overall Size2 :

2160mm (L) x 850mm (W) x 590mm (H).

Suitable size for mattress : 1905mm (L) x 810 mm (W)

#### OPTIONAL ACCESSORIES:

- Castors (Set of four) two with brake.
- Saline rod
- A pair of Side Railings.
- PU foam mattress.



### MW - 03(A): SEMI FOWLER BED (STANDARD)

Overall Size1 :

2060mm (L) x 925mm (W) x 590mm (H).

Suitable size for mattress : 1905mm (L) x 890mm (W)

#### Features:

- 60mm x 30mm CRCA tubular frame work.
- Two section perforated CRCA sheet top.
- Backrest operated by screw and lever mechanism with S.S. folding handle.
- 32mm Dia stainless steel bows.
- Four locations for I.V. Rod.
- Pretreated and epoxy powder coated.
- Size available as per users need.

Overall Size2 :

1985mm (L) x 850mm (W) x 590mm (H).

Suitable size for mattress : 1830mm (L) x 810mm (W)

#### OPTIONAL ACCESSORIES:

- Castors (Set of four) two with brake.
- Saline Rod.
- A pair of Side railings.
- Urine bag holder.
- Crib attachment.
- Colored metal bow panel.
- PU foam mattress.


**MW - 01: PLAIN HOSPITAL BED**

Overall Size1 :

2010mm (L) x 925mm (W) x 590mm (H).

Suitable size for mattress : 1905mm (L) x 890mm (W)

**Features:**

- 60mm x 30mm CRCA tubular frame work.
- 32mm dia stainless steel bows.
- Four locations for I.V. Rod.
- Pretreated and epoxy powder coated.
- Size available as per users need.

**OPTIONAL ACCESSORIES:**

- a) Castors Set of four (two with brake).
- b) Saline rod.
- c) A pair of side railings.
- d) Crib attachment.
- e) Urine bag holder.
- f) Coloured metal bow panel.


**MW - 02: HOSPITAL BED WITH MANUAL BACKREST**

Overall Size1 :

2010mm (L) x 925mm (W) x 590mm (H).

Suitable size for mattress : 1905mm (L) x 890mm (W)

**Features:**

- 60mm x 30mm CRCA tubular frame work.
- Two section perforated CRCA sheet top.
- Manually operated six position backrest on ratchet.
- 32mm dia stainless steel bows.
- Four locations for I.V. Rod.
- Pretreated and epoxy powder coated.
- Size available as per users need.



Overall Size2 :

1935mm (L) x 850mm (W) x 590mm (H).

Suitable size for mattress : 1830mm (L) x 810mm (W)

**OPTIONAL ACCESSORIES:**

- a) Castors Set of four (two with brake).
- b) Saline rod.
- c) A pair of side railings.
- d) Crib attachment.
- e) Urine bag holder.
- f) Coloured metal bow panel.

# PATIENT TRANSFER EQUIPMENTS

## MW - 19 (A) STRETCHER TROLLEY (FOLDING)

Overall approx Size: L 1830mm x W 570mm x H 800mm

Folding approx Size: L 1000mm x W 570mm x H 800mm

### Features:

- Both side folding stainless steel top fixed on CRCA tube frame.
- Base Frame: 50mm x 25mm CRCA Tube.
- Mounted on 150mm Dia heavy duty castors.
- Pretreated & Epoxy powder coated.
- Also available in Stainless Steel Top & Stainless Steel Base.



## MW - 20 STRETCHER ON TROLLEY

Overall approx Size:

L 2030mm x W 610mm x H 870mm

### Features:

- Base Frame: 50mm x 25mm CRCA Tube.
- Removable Stainless Steel top fixed on CRCA tube frame.
- Mounted on 150mm Dia heavy duty castors.
- Pretreated & Epoxy Powder Coated.
- Also available in Stainless Steel.



# PATIENT TRANSFER EQUIPMENTS

## MW - 21 (A) EMERGENCY & RECOVERY TROLLEY (CASUALTY TROLLEY) Manually Height Adjustable:



L 1980mm x W 710mm x Height adjustable from 700mm to 890mm

### Features:

- Two sections removable & X-Ray permeable stretcher top.
- X-Ray cassette holder.
- Height adjustment & Trendelenburg / Reverse Trendelenburg position by screw & lever mechanism.
- Oxygen cylinder cage.
- I. V. Rod.
- S.S. swing away type side railings.
- Four heavy duty swivel castors of 125mm Dia. (two with brakes).
- Buffers at four corners.
- Pretreated & Epoxy powder coated.

## MW - 21 (B) EMERGENCY & RECOVERY TROLLEY Manually Adjustable by Hydraulic Pump



L 1980mm x W 710mm x Height adjustable from 660mm to 910mm

### Features:

- Two sections removable & X-Ray permeable stretcher top.
- X-Ray cassette holder.
- Height adjustment by foot operated Hydraulic pump.
- Gas spring assisted Trendelenburg / Reverse Trendelenburg positions.
- Oxygen cylinder cage.
- I. V. Rod.
- S. S. swing away type side railings.
- Four heavy duty swivel castors of 125mm Dia. (two with brakes).
- Buffers at four corners.
- Pretreated & Epoxy powder coated.

**MW - 12 BED SIDE LOCKER (STANDARD)**

Overall Approx Size: W 370mm x D 305mm x H 760mm

**Features:**

- Two Stainless Steel tops with three side raised & front side down.
- Railings on three sides above main cabinet.
- Handle & magnetic latch on door.
- CRCA sheet box mounted on 20mm x 20mm CRCA tubes frames.
- Pretreated & Epoxy powder coated


**MW - 13 BED SIDE LOCKER (DELUXE)**

Overall Approx Size : W 395 mm x D 370 mm x H 805 mm

**Features:**

- CRCA Sheet main cabinet and drawer.
- Mounted on 50mm dia. castors.
- Stainless Steel top with three sides raised & front side down.
- Buffers at rear side.
- Pretreated & Epoxy powder coated


**VS - MI - 14 BED SIDE STAND FULL S.S.**

Overall Approx Size: W 390mm x D 360mm x H 770mm

**Features:**

- Three Stainless Steel Shelves.
- Mounted on Stainless Steel tubular frame.
- Optional: 50mm dia. castors.


**MW - 39 CRASH CART**

Overall approx Size : W 970mm x D 490mm x H 1630mm

**Features:**

- Six drawer modular system
- Upper Two drawer with medicine containers.
- Detachable plastic bins.
- Oxygen cylinder cage.
- Top tray for monitor & pulse oximeter.
- S.S. I.V. rod & cardiac massage board.
- System mounted on 100mm dia castors (two with brake).
- Available in Stainless Steel / Powder coated frame with S.S. top.


**MW - 50 ECG MACHINE TROLLEY**

Size : W 355mm x D 480mm x H 665mm

**Features:**

- 3 trays, 50mm dia wheels and Cable Hanger.
- Pre-treated & Epoxy powder coated.
- Available in Stainless Steel / Powder coated frame with S.S. top.


**MW - 17 HEIGHT ADJUSTABLE BED SIDE TABLE**

 Overall approx Size:  
 L 760mm x W 360mm

**Features:**

- Height Adjustable from 740mm to 1040mm
  - Stainless Steel top.
  - Height adjustable by bevel gear mechanism.
  - Mounted on 50mm dia castors.
  - Pre-treated & Epoxy powder coated.
- OPTIONAL** : Height Adjustable by gas spring.


**MW - 45 WALL MOUNTING MONITOR STAND**

Tray Size: W 270mm x D 380mm

- Pre-treated & epoxy powder coated.

**MW - 43 GYNAEC EXAMINATION COUCH**

Overall approx Size: L 1830mm x W 575mm x H 850mm

**Features:**

- Frame: CRCA Rectangular Tubes, CRCA sheet cabinets.
- Three cabinets one with lock & drawer with lock
- 100 mm thick reversible foam mattress with "U" cut.
- Examination Trolley with three trays.
- Sliding step stool & Writing pad.
- Pre-treated & Epoxy powder coated.
- OPTIONAL:** Lock facility for remaining Cabinets.


**MW - 44 (A) GYNAEC EXAMINATION TABLE**

Overall approx. Size : L 1820mm x W 560mm x H 815mm

**Features:**

- Frame: Press bend CRCA sheet Top with backrest adjustable on ratchet.
- Perineal "U" cut at foot end.
- CRCA tubular legs.
- Step stool & Mattress.
- Lithotomy bars with stirrups.
- Pre-treated & Epoxy powder coated

**MW - 44 (B) PLAIN EXAMINATION TABLE**

Overall approx Size: L 1820mm x W 560mm x H 815mm

**Features:**

- Frame: Press bend CRCA sheet Top with backrest adjustable on ratchet.
- CRCA tubular legs.
- Step stool & Mattress.
- Pre-treated & Epoxy powder coated





### MW - 42 (A) EXAMINATION COUCH

Overall Approx Size: L 1830mm x W 575mm x H 850mm

#### Features:

- Frame: CRCA Rectangular Tubes, CRCA sheet cabinets
- 100 mm thick reversible foam mattress.
- Provided with 4 cabinets, 1 with lock & a drawer with lock
- Sliding stepping stool & Writing pad
- Pre-treated & epoxy powder coated
- Optional: Lock facility for remaining Cabinets.



### MW - 42 (B) EXAMINATION COUCH WITH GAS SPRING OPERATED BACKREST

Overall Approx Size: L 1830mm x W 575mm x H 850mm

#### Features:

- Frame: CRCA Rectangular Tubes, CRCA sheet cabinets
- Sliding stepping stool.
- Backrest on Gas Spring
- Provided with Three drawers & three cabinets.
- One cabinet & One drawer with lock
- 100 mm thick reversible foam mattress
- Pre-treated & epoxy powder coated.
- Optional: Lock facility for remaining Drawers & Cabinets



### MW - 42 (C) SIMPLE EXAMINATION COUCH

Overall approx Size: L 1830mm x W 575mm x H 850mm

#### Features:

- Frame: CRCA Rectangular Tubes.
- CRCA sheet cabinet
- 100 mm thick reversible foam mattress
- One cabinet with lock, drawer with lock
- Sliding stepping stool & Writing pad
- Pre-treated & epoxy powder coated
- Accessories: Stepping stools & Mattress

**MW - 33 MONITOR TROLLEY**

Overall Size: W 965mm x D 660mm x H 1498mm  
Shelves size: W 610mm x D 510mm

**Features:**

- Adjustable shelves.
- Mounted on 125mm Dia. 125mm dia Heavy duty castors : 2 with brakes
- Eight point electrical board.
- Cylinder cage
- Available in Stainless Steel / Powder Coated


**MW - 38 INSTRUMENT CABINET**

Overall Size: W 610mm x D 355mm x H 1535mm

**Features:**

- Rectangular tubular structure.
- Sheet Panel on back side.
- Five shelves, door & two side of 5mm thick fine glass.
- Available in Stainless Steel / Powder Coated


**MW 23 C- OPERATION TABLE (ALL FUNCTION)**

Overall approx size: L 1830mm x W 710mm x Ht. Adjustable from 735mm to 970mm

**Features:**

- Stainless steel top and Stainless Steel Base
- Kidney bridge.
- Removable & folding foot end with drainage S.S. tray.
- Height adjustment and head low/up (T & RT) position by screw and lever mechanism.
- Lateral Tilt Position
- Lithotomy crutches.
- Gear type angular arrangement for lithotomy bars.
- Mattress, Step Stool & Arm boards.
- Stainless Steel removable side railings.

**MW - 28 INSTRUMENT TROLLEY**

(A) INSTRUMENT TROLLEY Overall approx.  
Size: L 710mm x B 460mm x H 850mm  
(B) INSTRUMENT TROLLEY Overall approx.  
Size: L 965mm x B 460mm x H 850 mm

**Features:**

- Stainless steel square tubes.
- Stainless Steel top & bottom shelves with 3 side railing on top shelf.
- Trolley mounted on 100mm Dia. castors.
- Available in full Stainless Steel / powder coated frame with Stainless Steel shelves.
- Also available in round tube frame. Supplied in KDC.


**MW - 29 MAYO'S TROLLEY (SINGLE BAR)**

Tray Size: L 555mm x W 420mm  
Height Adjustable From 845mm to 1300mm

**Features:**

- Frame: Stainless Steel Rectangular & Square Tubes.
- Trolley mounted on 50mm Dia. castors.
- Available in full Stainless Steel / powder coated frame with S.S. Tray.


**MW - 30 MAYO'S TROLLEY (DOUBLE BAR)**

Tray Size: L 770mm x W 510mm  
Height Adjustable From 845mm to 1300mm

**Features:**

- Frame: Stainless Steel Rectangular & Square Tubes.
- Trolley mounted on 50mm Dia. castors.
- Available in full Stainless Steel / powder coated frame with S.S. Tray.

**MW - 48 (A) WALL MOUNTED SCRUB**

Size : W 482mm x D 457mm

**Features:**

- Length as per requirement.
- Fully Stainless Steel.
- Provided with instrument washing tray.

OPTIONAL ACCESSORIES AT EXTRA COST:

- a) Stainless Steel sleeper stand.
- b) Long handle tap.

**MW - 48 (B) FOOT OPERATED SCRUB**

Size : W 558/mm x D 457mm x 1245mm

**Features:**

- Length as per requirement.
- Fully Stainless Steel.
- Instrument washing tray.




**A**

**B**
**MW - 18 DRUG TROLLEY**

 A) Size : W 490mm x D 580mm x H 1040 mm,  
4 trays & 100 mm dia castors

B) Size : W 355mm x D 480mm x H 665mm

**Features:**

- 3 trays & 50mm Dia. castors.
- Preatreated & Epoxy powder coated.
- Available in full Stainless Steel / Epoxy powder coated.
- Both are also available in Stainless Steel.

**MW - 40 (A) BABY RESUSCITATION TROLLEY**

Overall approx Size: W 390mm x D 800mm x H 1080mm

**Features:**

- Stainless Steel tray in head position.
- Three drawer modular system.
- Adjustable I.V. Rod.
- Oxygen cylinder cage.
- Adjustable Examination Lamp.
- Bottom tray for mounting suction apparatus.
- Trolley mounted on 100 mm Dia castors.
- Preatreated & Epoxy powder coated.
- Available in full Stainless Steel / powder coated.


**MW - 37 DRESSING TROLLEY WITH BOWL & BUCKET**

Overall trolley Size: L 1080 x B 460mm x H 850mm

**Features:**

- Frame: Stainless Steel square tubes of 25mm x 25mm.
- Two Stainless Steel shelves with railing top.
- Stainless Steel bowl Dia 345mm, bucket Dia 360mm.
- Available in full Stainless Steel / Powder coated frame with Stainless Steel shelves & railings.


**MW - 31 REVOLVING STOOL**

Overall approx Size : Top 300mm Dia with  
Height Adjustable 465mm to 710mm Dia with

**Features:**

- (A) Stainless Steel top & Epoxy powder coated base.
  - (B) Stainless Steel Top & Stainless Steel base
  - (C) Stainless Steel Top & Stainless Steel base with Backrest.
- Optional : • Height adjustment by Gas spring • Cushions for seat & Backrest.


**MW - 34 SOILED LINEN TROLLEY**

Linen Bag Size: 510mm dia x H 780mm

**Features:**

- Trolley mounted on 100mm Dia. three castors.
- Available in full Stainless Steel / Epoxy powder coated.


**MW - 41 OXYGEN CYLINDER TROLLEY**

- A) Pull Type with three castors 50mm Dia.
- B) Push Type with two castors 100mm Dia.
- C) Jumbo Oxygen Trolley with two castors 125mm Dia.


**MW - 36 (A) TWO TIER BOWL STAND**

Overall approx Size: W 405mm x D 405mm x H 840mm

**Features:**

- Two Stainless Steel bowls dia 360mm each.
- Available in Stainless Steel.
- Epoxy powder coated.


**MW - 35 KICK BUCKET**

Overall approx Size: W 365mm x D 365mm x H 330mm

**Features:**

- Two Stainless Steel bowls dia 360mm each.
- Available in Stainless Steel.
- Epoxy powder coated.


**MW - 11 (A) CRIB ATTACHMENT**

**MW - 11 (B) CRIB WITH STAND**

**MW - 17 (C) FOOD TABLE-FULL S.S.**  
 Stainless Steel frame and Stainless Steel top.  
 Overall approx  
 Size: L 680 mm x W 360 mm x H 310 mm

**MW - 10 (A) I.V. ROD FIX HEIGHT**  
 Stainless Steel

**MW - 10 (B) I.V. ROD TELESCOPIC**  
 Stainless Steel

**MW - 32 SALINE STAND**  
 Ht. Adjustable from Min 1370 mm Max : 2435 mm  
 • Available in Full Stainless Steel / Epoxy powder coated stand & Stainless Steel saline Rod.

**MW - 16 BED SIDE SCREEN**  
 Overall approx. Size L open 2360mm / closed 1200mm x W 550mm x H 1730mm  
 • CRCA tubular structure mounted on 50mm dia wheels two sides folding.

**MW - 15 MULTIPURPOSE STOOL**  
 Overall Top 280mm x 280mm & H 445mm  
 Frame: 25mm x 25mm, CRCA Square Tube  
 • Stainless Steel Top & Epoxy powder coated Base.  
 • Also Available in Stainless Steel Top & Stainless Steel base.

**MW - 46 BED ELEVATOR**


**MW - 49 RECLINING CHAIR**

Overall Approx Size: 1730mm (L) x 558mm (W)  
 Chair Position: 780mm (L) x 558mm (W)

- Sliding with folding mattress
- Two locations fo


**MW- 09 (A): ATTENDANT BED**

Pretreated and Epoxy powder coated.  
 Overall Size1 :  
 1830mm (L) x 610mm (W) x 470mm (H).  
 50mm x 25mm CRCA rectangular tubular frame work.

**MW 06 PEADIATRIC BED WITH DROP SIDE RAILINGS**

Overall Size1 :  
 L 1450mm x W 840mm x H 590mm  
 Mattress size: L 1350mm x W 750mm

- Two locations for I.V. Rod.
- Perforated CRCA Sheet Top
- Full length drop side rails.
- Pretreated & Epoxy powder coated.



## ELECTRO HYDRAULIC O.T. TABLE WITH LONGITUDINAL TOP SLIDE & MANUAL OVERRIDE



**MW-TYCO 8006**

**TECHNICAL DATA**

Top Slide both Head & leg side : 30 cm manual operated  
 Minimum Height : 75 cm  
 Maximum Height : 110 cm  
 Length of Table Top : 215 cm  
 Width of Table Top : 50 cm  
 Lateral Tilt : + 20° / -20°  
 Trendelenburg : +25° / -25°  
 Back Section: +90° / -40° Head & leg Section: +90° / -90° manual operated  
 Kidney Elevator : 15 cm manual operated

**SALIENT FEATURES:**

100% Life Saving due to low voltage Electricity Supply.  
 Offline Battery Back-Up.  
 800 lb. weight capacity.  
 Full Table Functions, Full Articulation Capability.  
 Full Body imaging without special accessories.  
 Durable one-touch pendant control.  
 Casters lockable with single action brake  
 All stainless steel accessories  
 Electrically Controlled Height, Trendelenburg/Reverse, Lateral Tilt.  
 Radiolucent table top  
 Support on a heavy base column covered with stainless steel sheet.

**ACCESSORIES OF O.T. TABLE**

Anaesthetic Frame (1Pcs)  
 Padded Side Supports (1 Pair)  
 Intravenous Armboard (1 Pair)  
 Padded Shoulder Supports (1 Pair)  
 Knee Crutches (1 Pair)  
 One Matters Set Thickness 50mm

### Various positions of MW-TYCO 8006



## ELECTRO HYDRAULIC O.T. TABLE WITH LONGITUDINAL & LATERAL TOP SLIDE

**MW-TYCO 8007**



**TECHNICAL DATA**

Longitudinal Top Slide : 40 cm  
 Lateral Top Slide : 20 cm  
 Minimum Height : 75 cm  
 Maximum Height : 100 cm  
 Length of Table Top : 215 cm  
 Width of Table Top : 50 cm  
 Lateral Tilt : + 20° / -20°  
 Trendelenburg : +25° / -25°

**Specially Design  
for  
Vascular Surgery**

**SALIENT FEATURES:**

- 100% Life Saving due to low voltage Electricity Supply.
- Offline Battery Back-Up.
- 800 lb. weight capacity.
- Full Table Functions, Full Articulation Capability.
- Full Body imaging without special accessories.
- Durable one-touch pendant control.
- Casters lockable with single action brake
- All stainless steel accessories
- Electrically Controlled Height, Trendelenberg/Reverse, Lateral Tilt.
- Radiolucent table top
- Support on a heavy base column covered with stainless steel sheet.

**ACCESSORIES OF O.T. TABLE**

- Anaesthetic Frame (1Pcs)
- Padded Side Supports (1 Pair)
- Intravenous Armboard (1 Pair)
- Padded Shoulder Supports (1 Pair)
- Knee Crutches (1 Pair)
- One Matters Set Thickness 50mm

**Various positions of MW-TYCO 8007**



## C - ARM PROVISION HYDRAULIC O.T. TABLE WITH LONGITUDINAL TOP SLIDE



**MW-TYCO 8004**

### TECHNICAL DATA

Top Slide : 20 cm  
 Minimum Height : 80 cm  
 Maximum Height : 110 cm  
 Length Table Top : 195 cm to 215 cm  
 Width of Table Top : 48 cm  
 Lateral tilt :  $20^{\circ} \pm 2^{\circ}$   
 Trendelenburg :  $30^{\circ} \pm 2^{\circ}$   
 Chair Position :  $-22^{\circ} / +90^{\circ}$   
 Head Section:  $+90^{\circ} / -90^{\circ}$   
 Leg Position :  $+90^{\circ} / -45^{\circ}$   
 Kidney Elevation: 15 cm

### SALIENT FEATURES:

Compatible with 'C'-Arm Image Intensifier used orthopaedics surgery.  
 Can also be used for general surgeries like Cardiac, Urology, Endoscopy.  
 The table top is made of radio-translucent material.  
 Head and leg sections are Interchangeable.  
 Kidney bridge detachable or inbuilt as per customer requirement.  
 Base and center column are made of stainless steel or Fiber Glass for easy cleaning and hygienic.  
 Table framed, Side rails and Accessories are made up of stainless steel.

### ACCESSORIES OF O.T. TABLE

Anaesthetic Frame (1Pcs)  
 Padded Side Supports (1 Pair)  
 Intravenous Armboard (1 Pair)  
 Padded Shoulder Supports (1 Pair)  
 Knee Crutches (1 Pair)  
 One Matters Set Thickness 50mm

### Various positions of MW-TYCO 8004



## C - ARM PROVISION HYDRAULIC O.T. TABLE



**MW-TYCO 8004A**

### TECHNICAL DATA

Minimum Height : 75 cm  
 Maximum Height : 105 cm  
 Length Table Top : 195 cm to 215 cm  
 Width of Table Top : 48 cm  
 Lateral tilt :  $20^{\circ} \pm 2^{\circ}$   
 Trendelenburg :  $30^{\circ} \pm 2^{\circ}$   
 Chair Position :  $-22^{\circ} / +90^{\circ}$   
 Head Section :  $+90^{\circ} / -90^{\circ}$   
 Leg Position :  $+90^{\circ} / -45^{\circ}$   
 Kidney Elevation: 15 cm

### SALIENT FEATURES:

Compatible with 'C'-Arm Image Intensifier used orthopaedics surgery.  
 Can also be used for general surgeries like Cardiac, Urology, Endoscopy.  
 The table top is made of radio-translucent material.  
 Head and leg sections are Interchangeable.  
 Kidney bridge detachable or inbuilt as per customer requirement.  
 Base and center column are made of stainless steel or Fiber Glass  
 for easy cleaning and hygienic.  
 Table framed, Side rails and Accessories are made up of stainless steel.

### ACCESSORIES OF O.T. TABLE

Anaesthetic Frame (1Pcs)  
 Padded Side Supports (1 Pair)  
 Intravenous Armboard (1 Pair)  
 Padded Shoulder Supports (1 Pair)  
 Knee Crutches (1 Pair)  
 One Matters Set Thickness 50mm

### Various positions of MW-TYCO 8004A



# MAJOR HEAD END CONTROL HYDRAULIC O.T. TABLE AS PER ISI - 5291



**MW-TYCO 8003**

**TECHNICAL DATA**

Minimum Height :  $75 \pm 7.5$  cm  
 Maximum Height :  $115 \pm 5$   
 Length Table Top : 180 cm to 195 cm  
 Width of Table Top : 50 cm  
 360° Revolving  
 Lateral tilt :  $20^\circ \pm 2^\circ$   
 Trendelenburg :  $45^\circ \pm 3^\circ$   
 Chair Position :  $-22^\circ / +90^\circ$   
 Flex/Reflex :  $90^\circ / 22^\circ$   
 Kidney Position : 15 cm

**SALIENT FEATURES:**

Versatile Operation table for General Surgical Procedures.  
 Accurate positions can be achieved with easy mechanical operation  
 Five section Stainless Steel top  
 Table top can be rotate 360° on central axis  
 Heavy duty Single Cylindrical hydraulic pump  
 Base and Column covered with stainless steel fittings for hygiene and durability  
 Movable on High Quality Castors  
 Equipped with base locking mechanism for firm and rigid locking of the table detachable leg section

**ACCESSORIES OF O.T. TABLE**

Anaesthetic Frame (1Pcs)  
 Padded Side Supports (1 Pair)  
 Intravenous Armboard (1 Pair)  
 Padded Shoulder Supports (1 Pair)  
 Knee Crutches (1 Pair)  
 One Matters Set Thickness 50mm  
 Body Strap

**Various positions of MW-TYCO 8003**



# MINOR HEAD END CONTROL HYDRAULIC O.T. TABLE



**MW-TYCO 8002**

**TECHNICAL DATA**

Minimum Height :  $75 \pm 7.5$  cm  
 Maximum Height :  $115 \pm 5$   
 Length Table Top : 180 cm to 195 cm  
 Width of Table Top : 50 cm  
 Lateral tilt :  $20^\circ \pm 2^\circ$   
 Trendelenburg :  $45^\circ \pm 3^\circ$   
 Chair Position :  $-22^\circ / +90^\circ$   
 Flex/Reflex :  $90^\circ / 22^\circ$   
 Kidney Position : 15 cm


**SALIENT FEATURES:**

Single Action Hydraulic Lift  
 Heavy Mild Steel Base with Fibre wheels provides stability to the entire table  
 Auto floor locking device  
 The detachable leg section provides easy lithotomy Positioning  
 The improved Hydraulic system provides smooth movement and the speed of lowering the table is permanently set.

**ACCESSORIES OF O.T. TABLE**

Padded Side Supports (1 Pair)  
 Intravenous Armboard (1 Pair)  
 Padded Shoulder Supports (1 Pair)  
 Knee Crutches (1 Pair)  
 One Matters Set Thickness 50mm

**Various positions of MW-TYCO 8002**



## SIDE END CONTROL HYDRAULIC O.T. TABLE



**MW-TYCO 8001**

**TECHNICAL DATA**

Minimum Height : 75 cm  
 Maximum Height : 95 cm  
 Length Table Top : 210 cm  
 Width of Table Top : 50 cm  
 Lateral tilt :  $20^{\circ} \pm 2^{\circ}$   
 Trendelenburg :  $30^{\circ} \pm 3^{\circ}$   
 Hydraulic Lift: 22 cm  
 Back Position:  $+90^{\circ} / -40^{\circ}$   
 Kidney Position : 12 cm

**SALIENT FEATURES:**

Single Action Hydraulic Lift  
 Heavy Mild Steel Base with Fibre wheels provides stability to the entire table  
 Auto floor locking device  
 The detachable leg section provides easy lithotomy Positioning  
 The improved Hydraulic system provides smooth movement and the speed of lowering the table is permanently set.

**ACCESSORIES OF O.T. TABLE**

Padded Side Supports (1 Pair)  
 Intravenous Armboard (1 Pair)  
 Padded Shoulder Supports (1 Pair)  
 Knee Crutches (1 Pair)  
 One Matters Set Thickness 50mm

**Various positions of MW-TYCO 8001**



## 360° HYDRAULIC O.T. TABLE as per ISI - 6106



### TECHNICAL DATA

Minimum Height :  $75 \pm 5$  cm  
 Maximum Height :  $95 \pm 5$  cm  
 Length Table Top : 180 cm to 195 cm  
 Width of Table Top : 50 cm  
 Lateral tilt :  $20^\circ + 2^\circ$   
 Trendelenburg :  $30^\circ + 3^\circ$   
 Hydraulic Lift: 22 cm


### SALIENT FEATURES:

Single Action Hydraulic Lift  
 Heavy Mild Steel Base with Fibre wheels provides stability to the entire table  
 Auto floor locking device  
 The detachable leg section provides easy lithotomy Positioning  
 The improved Hydraulic system provides smooth movement and the speed of lowering the table is permanently set.


### ACCESSORIES OF O.T. TABLE

Anaesthetic Frame (1Pcs)  
 Padded Side Supports (1 Pair)  
 Intravenous Armboard (1 Pair)  
 Padded Shoulder Supports (1 Pair)  
 Knee Crutches (1 Pair)  
 One Matters Set Thickness 50mm

### Various positions of MW-TYCO 8001R



## REMOTE CONTROL ELECTRO O.T. TABLE WITH LONGITUDINAL TOP SLIDE



**MW-TYCO 8005**

**TECHNICAL DATA**

Top Slide both Head & leg side : 30 cm  
 Minimum Height : 75 cm  
 Maximum Height : 105 cm  
 Length of Table Top : 215 cm  
 Width of Table Top : 50 cm  
 Lateral Tilt : + 20° / -20°  
 Trendelenburg : +25° / -25°  
 Back Section: +90° / -40°  
 Head & leg Section: +90° / -90° manual operated  
 Kidney Elevator : 15 cm manual operated


**SALIENT FEATURES:**

100% Life Saving due to low voltage Electricity Supply.  
 Offline Battery Back-Up.  
 800 lb. weight capacity.  
 Full Table Functions, Full Articulation Capability.  
 Removable Leg section.  
 Full Body imaging without special accessories.  
 Durable one-touch pendant control.  
 Electrically Controlled Height, Trendelenburg/Reverse, Lateral Tilt.  
 4-Section seamless moulded mattress  
 Support on a heavy base column covered with stainless steel sheet.

**ACCESSORIES OF O.T. TABLE**

Anaesthetic Frame (1Pcs)  
 Padded Side Supports (1 Pair)  
 Intravenous Armboard (1 Pair)  
 Padded Shoulder Supports (1 Pair)  
 Knee Crutches (1 Pair)  
 One Matters Set Thickness 50mm

### Various positions of MW-TYCO 8005



**(MODEL - C.S. 19D+19D)**  
**CEILING SHADOWLESS O.T. LIGHT**  
**TWIN DOME WITH DIACHROIC REFLECTOR**



## TECHNICAL SPECIFICATION

DOME DIA : 22" + 22"


HALOGEN BULB : 24V/150W

COLOR : 4200 K

POWER : 220V AC/DC

INTENSITY : 1,10,000 Lux + 1,10,000 Lux

PROVIDED WITH AUXILIARY BULB TO AVOID LIGHT FAILURE DURING  
SURGERY



**Imported Heat Observing UV Filter Glass**

**(MODEL - C.S. - 7+4)**  
**CEILING SHADOWLESS TWIN O.T. LIGHT**  
**SEVEN & FOUR REFLECTOR**



## TECHNICAL SPECIFICATION

DOME DIA : 70CM + 55CM

HALOGEN BULB : 12V - 55W X 11


ACTION RADIUS : 160 CM

FIELD AREA : 20CM

COLOR : 4200 K

POWER : 220V AC/DC

INTENSITY : 1,00,000 Lux + 70,000 Lux



**Imported Heat Observing UV Filter Glass**

**(MODEL - C.S. - 4+4)**  
**CEILING SHADOWLESS O.T. LIGHT**  
**FOUR REFLECTOR TWIN**



## TECHNICAL SPECIFICATION

DOME DIA : 55CM + 55CM

HALOGEN BULB : 12V - 55W X 8


ACTION RADIUS : 160 CM

FIELD AREA : 20CM

COLOR : 4200 K

POWER : 220V AC/DC

INTENSITY : 70,000 Lux + 70,000 Lux



**Imported Heat Observing  
UV Filter Glass**

**(MODEL - C.S. -19 + 19)**  
**CEILING SHADOWLESS O.T. LIGHT**  
**TWIN DOME**



## TECHNICAL SPECIFICATION

DOME DIA : 48CM + 48CM

HALOGEN BULB : 24V/150W

COLOR : 4200 K

POWER : 220V AC/DC

INTENSITY : 50,000 Lux + 50,000 Lux



**Imported Heat Observing UV Filter Glass**

PROVIDED WITH AUXILIARY BULB TO AVOID LIGHT FAILURE DURING  
SURGERY

**(MODEL - C.S. 19D)**  
**CEILING SHADOWLESS O.T. LIGHT**  
**SINGLE DOME WITH DIACHROIC REFLECTOR**



## **TECHNICAL SPECIFICATION**

DOME DIA : 22"

HALOGEN BULB : 24V/150W

COLOR : 4200 K

POWER : 220V AC/DC

INTENSITY : 1,10,000 Lux

PROVIDED WITH AUXILIARY BULB TO AVOID LIGHT FAILURE DURING  
SURGERY




**Imported Heat Observing UV Filter Glass**

**(MODEL - C.S. - 7)**  
**CEILING SHADOWLESS O.T. LIGHT**  
**SINGLE DOME SEVEN REFLECTOR**



## **TECHNICAL SPECIFICATION**

DOMES DIA : 70CM  
HALOGEN BULB : 12V - 55W X 7  
ACTION RADIUS : 160 CM  
FIELD AREA : 20CM  
COLOR : 4200 K  
POWER : 220V AC/DC  
INTENSITY : 1,00,000 Lux



**Imported Heat Observing  
UV Filter Glass**

**(MODEL - C.S. - 4)**  
**CEILING SHADOWLESS O.T. LIGHT**  
**SINGLE DOME FOUR REFLECTOR**



## TECHNICAL SPECIFICATION

DOME DIA : 55CM

HALOGEN BULB : 12V - 55W X 4


ACTION RADIUS : 160 CM

FIELD AREA : 20CM

COLOR : 4200 K

POWER : 220V AC/DC

INTENSITY : 70,000 Lux



**Imported Heat Observing  
UV Filter Glass**

**(MODEL - C.S. -19)**  
**CEILING SHADOWLESS O.T. LIGHT**  
**SINGLE DOME**



## **TECHNICAL SPECIFICATION**

DOME DIA : 48CM

HALOGEN BULB : 24V/150W

COLOR : 4200 K

POWER : 220V AC/DC

INTENSITY : 50,000 Lux

PROVIDED WITH AUXILIARY BULB TO AVOID LIGHT FAILURE DURING  
SURGERY



**Imported Heat Observing UV Filter Glass**

**(MODEL - M.S. 19D)**  
**MOBILE SHADOWLESS O.T. LIGHT**  
**SINGLE DOME WITH DIACHROIC REFLECTOR**



**Imported Heat Observing  
UV Filter Glass**



**TECHNICAL SPECIFICATION**

DOME DIA : 22"


HALOGEN BULB : 24V/150W

COLOR : 4200 K

POWER : 220V AC/DC

INTENSITY : 1,10,000 Lux

PROVIDED WITH AUXILIARY BULB TO AVOID LIGHT FAILURE DURING  
SURGERY

**(MODEL - M.S. - 4)****MOBILE SHADOWLESS O.T. LIGHT  
SINGLE DOME FOUR REFLECTOR**

Imported Heat Observing  
UV Filter Glass

**TECHNICAL SPECIFICATION**

DOME DIA : 55CM

HALOGEN BULB : 12V - 55W X 4

FIELD AREA : 20CM

COLOR : 4200 K

POWER : 220V AC/DC

INTENSITY : 70,000 Lux



**(MODEL - M.S. -19)**  
**MOBILE SHADOWLESS O.T. LIGHT**  
**SINGLE DOME**



**Imported Heat Observing  
UV Filter Glass**

**TECHNICAL SPECIFICATION**

DOME DIA : 48CM

HALOGEN BULB : 24V/150W

COLOR : 4200 K

POWER : 220V AC/DC

INTENSITY : 50,000 Lux

PROVIDED WITH AUXILIARY BULB TO AVOID LIGHT FAILURE DURING  
SURGERY



**(MODEL - M.S. 19DB)**  
**MOBILE SHADOWLESS O.T. LIGHT**  
**SINGLE DOME WITH DIACHROIC REFLECTOR**  
**AND BATTERY BACK-UP**



**Imported Heat Observing  
UV Filter Glass**

**TECHNICAL SPECIFICATION**

DOME DIA : 22"

HALOGEN BULB : 24V/150W

COLOR : 4200 K


POWER : 220V AC/DC

INTENSITY : 1,10,000 Lux

PROVIDED WITH AUXILIARY BULB TO AVOID LIGHT FAILURE DURING  
SURGERY



**(MODEL - M.S. - 4)**  
**MOBILE SHADOWLESS O.T. LIGHT**  
**SINGLE DOME FOUR REFLECTOR**  
**WITH BATTERY BACK-UP**



**Imported Heat Observing  
UV Filter Glass**

**TECHNICAL SPECIFICATION**

DOME DIA : 55CM

HALOGEN BULB : 12V - 55W X 4

FIELD AREA : 20CM

COLOR : 4200 K

POWER : 220V AC/DC

INTENSITY : 70,000 Lux





We introduce ourselves as one of the leading concern dealing in manufacturing and installing of lab furniture and clean room furniture, Fume Hood, along with office designing and furnishing. We are one of the finest laboratory furniture suppliers in india, We are working on C.N.C. Machines for METAL AND WOOD.

### HPLC TABLE IN INSTRUMENT LAB

Elegant, white coloured table top. All the electric trunking is suitably placed to take power. The height of Instrument benches is 750mm (sitting ht) and leg space is provided in between. Keyboard drawer and CPU shelves are adding the facility to work with computerized instruments.



### CORNER SINK & THREE WAY WATER TAP

Sink Benches are suitable for Laboratory consists of assembly of combination of storage modules with shutters with shelf inside and Tops of selected material. The main feature is to facilitate the Sink along with the water tap integrated with the work bench.

Three way water tap is made up of forged brass body with powder coated. The material of sink is non reactive to acids/alkalies/solvent/hot & cold water. The sink is easy to maintain and is rust proof. The table top is also easy to maintain.

### THE TECHNICAL SPECIFICATION WILL BE AS FOLLOWS.

1. All the metal used will be 1.00 mm thick G.I. STEEL, (The make will be TATA OR JINDAL)
2. 1.2mm G.I. Sheets will be used for Fume Hood where it is required.
3. 65-85 thick micron epoxy powder coating will be done for all the furniture. (The powder will be of Asian or Nerolac make)
4. Door will be of double walls. The inner wall and back of the storage along with bottom of the drawers. will be of composite material, which is not affected by any acid/alkalies/solvent/hot water/cold water and is rust proof.
5. All the sliding channels will be of telescopic channels of good quality.
6. All the Hinges which are introduced for the first time in India, does not require any maintenance. (DEVELOPED INHOUSE)
7. We will use Silica coated sink, which is non reactive to acids/alkaline/solvent/hot & cold water, & easy to maintain.
8. Table Top will be 18-/+ 16mm thick Jet Black Granite. (Approximate cost will be Rs.155/Sft)
9. All the sinks will be fitted from below the granite, as it will be easy to clean.
10. The working Height of the table in chemical lab will be 900mm, which can be changed up to 850 with adjustable screw, fixed below the frame of the table.
11. The working Height of the table in Instrument lab will be 750mm, which can be changed up to 700 with adjustable screw, fixed below the frame of the table.
12. The Fume Hoods will have automatic sash opening/closing facility.
13. 95% of the material which will be used in the Fume Hood will be non reactive to acids/Alkaline/solvent/hot and cold water and rust proof, and will be white in color and easy to maintain.
14. There will be no Hinges in storage area for the shutter for the Fume Hood.
15. Color combination VIP blue drawer & broken white body. Or as per party requirement.
16. Delivery period will be around four to six weeks from the date of receipt of Purchase order in our favor.



### INSTRUMENT WORKING TABLE

Modern lab with elegant table top. The table has all the requirements, such as elegant look, key board, electrical points, ready to use for any kind of instrument in any lab.

### ISLAND TABLE WITH ULTRA CLEAN GLASS TOP

Island Benches are suitable for chemical testing in a Laboratory. They consist of an assembly of combination of storage modules with drawers & shutters or only shutter with shelf inside and tops of selected material.

Open shelves are provided between each module which are very handy for keeping day to day samples, or keeping pipets, beakers etc. that can be removed immediately.



### REAGENT RACKS WITH GLASS TABLE TOP

They are fitted in the middle of island table operational both sides as well as working table fitted along the wall. The racks will have 2 shelves to keep the reagent made up of special white material easy to clean. The reagent rack is made up of special pillars with removable plates for electrical switches to be fitted on the plate as it is placed in corner to avoid explosive acid falling on the switch.

## ADVANTAGES OF OUR FUME HOOD

- Our is the only fume hood with automatic door operation with a click of a button.
- **FREE FROM BOLTS & CAPS** our fume hood are free from bolts & cap system, from inside & from outside also, so no question of rusting of nuts & bolts.
- **CRAFTSMANSHIP** excellent, since we are using only C. N. C. Machines & better than any of fume hood manufactured in India.
- **DESIGN** our fume hood are compatible with all the international standards may be better, hence no need to waste extra money in the name of ashare or any other standard.
- **INTERLINING** special chemical/ acid/ alkalis/fire/ & scratch resistance material for main chamber and for the storage cabinet.
- **ELECTRIC FITTINGS** in house almost fire proof type electrical tube light along with is make electron wirings
- Electrically operated door (sash). We are the only manufacturer who provide with automatic door operation without charging our customers extra.
- **NOZZLES & GAS VALVES** are color coded as per din 12920
- **HINGES** we are the only manufacturer in India who have eliminated hinges in the fume hood hence no question of rusting
- **SILENT BLOWER'S** smooth & silent blowers, highly efficient, that do not produce any noise & vibration.
- Easy to assemble and dismantle and to services.
- Level adjusting legs for + 10mm height adjustment.
- Shortest delivery period, best value for your money.