



“Brands of Genius are Admired,  
Brands of Wealth are Envied,  
Brands of Controversies are Feared;  
But only brands of Reputed Company are  
**TRUSTED.**”

**ACCUREX<sup>®</sup>**

**sero<sup>o</sup>react<sup>®</sup>**

# SeroReact Product Range

## SeroReact ASO/ RF/ CRP

- Latex agglutination slide test
- Qualitative and Semi-Quantitative testing possible
- Rapid result in 2 minutes
- Complete composite kit with controls
- Pack - Size : 25 Tests
- Analytical sensitivity optimized to:  
ASO - 200 IU/ml  
RF - 8 IU/ml  
CRP - 0.6 mg/dl

## SeroReact RPR

- Carbon flocculation slide test
- Qualitative and Semi-Quantitative testing possible
- Rapid result in 4 minutes
- Complete composite kit with controls
- Pack - Size : 50 Tests

## SeroReact WIDAL

- Qualitative and Quantitative slide test
- Detection of 'O', 'H', A(H) & B(H)
- Rapid result in 1 minute
- User-friendly color-coded packaging
- Pack - Size : 4 x 5 ml - for detection of 'O', 'H', A(H) & B(H)  
(2+2) x 5 ml - for detection of 'O' & 'H'

SeroReact

Assurance of Persistent Quality  
from **ACCUREX**  
A trusted name in Diagnostics

