



THE NEW BENCHMARK IN HEALTHCARE

ELECTRO HYDRAULIC O.T. TABLE WITH LONGITUDINAL & LATERAL TOP SLIDE



TYCO 8007

Specially Design
for
Vascular Surgery

TECHNICAL DATA : Longitudinal Top Slide : 40 cm
Lateral Top Slide : 20 cm
Minimum Height : 75 cm
Maximum Height : 105 cm
Length of Table Top : 215 cm
Width of Table Top : 50 cm
Lateral Tilt : + 20° / -20°
Trendelenburg : +25° / -25°

SALIENT FEATURES : 100% Life Saving due to low voltage Electricity Supply.
Offline Battery Back-Up.
800 lb. weight capacity.
Full Table Functions, Full Articulation Capability.
Full Body imaging without special accessories.
Durable one-touch pendant control.
Casters lockable with single action brake
All stainless steel accessories
Electrically Controlled Height , Trendelenburg/Reverse, Lateral Tilt.
Radiolucent table top
Support on a heavy base column covered with stainless steel sheet.

ACCESSORIES OF O.T. TABLE : Anaesthetic Frame (1Pcs)
Padded Side Supports (1 Pair)
Intravenous Armboard (1 Pair)
Padded Shoulder Supports (1 Pair)
Knee Crutches (1 Pair)
One Matters Set Thickness 50mm



THE NEW BENCHMARK IN HEALTHCARE

VARIOUS POSITIONS OF TYCO 8007



Trendelenburg



Lateral Top Slide



Reverse Trendelenburg



Lateral Tilt (Both)



Lateral Tilt Left



Lateral Tilt Right



Longitudinal Top Slide

ELECTRO HYDRAULIC O.T. TABLE WITH LONGITUDINAL TOP SLIDE & MANUAL OVERRIDE



TYCO 8006

TECHNICAL DATA : Top Slide both Head & leg side : 30 cm manual operated
Minimum Height : 75 cm
Maximum Height : 110 cm
Length of Table Top : 215 cm
Width of Table Top : 50 cm
Lateral Tilt : + 20° / -20°
Trendelenburg : +25° / -25°
Back Section: +90° / -40°
Head & leg Section: +90° / -90° manual operated
Kidney Elevator : 15 cm manual operated

SALIENT FEATURES : 100% Life Saving due to low voltage Electricity Supply.
Battery Back-Up.
800 lb. weight capacity.
Full Table Functions, Full Articulation Capability.
Removable Leg section.
Full Body imaging without special accessories.
Durable one-touch pendant control.
Casters lockable with single action brake
All stainless steel accessories
Electrically Controlled Height , Backrest, Trendelenburg/Reverse, Lateral Tilt.
Radiolucent table top
Manual over ride by rotating crank on base.
Head, Backrest, Trunk and Leg 4-section tabletop
5-Section seamless moulded mattress
Support on a heavy base column covered with stainless steel sheet.

ACCESSORIES OF O.T. TABLE : Anaesthetic Frame (1Pcs)
Padded Side Supports (1 Pair)
Intravenous Armboard (1 Pair)
Padded Shoulder Supports (1 Pair)
Knee Crutches (1 Pair)
One Matters Set Thickness 50mm



THE NEW BENCHMARK IN HEALTHCARE

VARIOUS POSITIONS OF TYCO 8006



Trendelenburg



Lateral Tilt (Both)



Reverse Trendelenburg



Chair Position



Flex Position



Reflex Position



Manual Selector for position



Kidney Position



Manual Lever for positions

REMOTE CONTROL ELECTRO O.T. TABLE WITH LONGITUDINAL TOP SLIDE



TYCO 8005

TECHNICAL DATA :

- Top Slide both Head & leg side : 30 cm
- Minimum Height : 75 cm
- Maximum Height : 105 cm
- Length of Table Top : 215 cm
- Width of Table Top : 50 cm
- Lateral Tilt : + 20° / -20°
- Trendelenburg : +25° / -25°
- Back Section: +90° / -40°
- Head & leg Section: +90° / -90° manual operated
- Kidney Elevator : 15 cm manual operated

SALIENT FEATURES :

- 100% Life Saving due to low voltage Electricity Supply.
- Offline Battery Back-Up.
- 800 lb. weight capacity.
- Full Table Functions, Full Articulation Capability.
- Removable Leg section.
- Full Body imaging without special accessories.
- Durable one-touch pendant control.
- Electrically Controlled Height , Backrest, Trendelenburg/Reverse, Lateral Tilt.
- Radiolucent table top
- 4-Section seamless moulded mattress
- Support on a heavy base column covered with stainless steel sheet.

ACCESSORIES OF O.T. TABLE :

- Anaesthetic Frame (1Pcs)
- Padded Side Supports (1 Pair)
- Intravenous Armboard (1 Pair)
- Padded Shoulder Supports (1 Pair)
- Knee Crutches (1 Pair)
- One Matters Set Thickness 50mm




THE NEW BENCHMARK IN HEALTHCARE

VARIOUS POSITIONS OF TYCO 8005



Trendelenburg



Lateral Tilt (Both)



Reverse Trendelenburg



Chair Position



Flex Position



Reflex Position



Kidney Position

C - ARM PROVISION HYDRAULIC O.T. TABLE



TYCO 8004 A

TECHNICAL DATA :

- Minimum Height : 75 cm
- Maximum Height : 105 cm
- Length Table Top : 195 cm to 215 cm
- Width of Table Top : 48 cm
- Lateral tilt : $20^{\circ} \pm 2^{\circ}$
- Trendelenburg : $30^{\circ} \pm 2^{\circ}$
- Chair Position : $-22^{\circ} / +90^{\circ}$
- Head Section : $+90^{\circ} / -90^{\circ}$
- Leg Position : $+90^{\circ} / -45^{\circ}$
- Kidney Elevation : 15 cm

SALIENT FEATURES :

- Compatible with 'C'-Arm Image Intensifier used for orthopaedics surgery. Can also be used for general surgeries like Cardiac, Urology, Endoscopy. The table top is made of radio-translucent material.
- Head and leg sections are interchangeable.
- Kidney bridge detachable or inbuilt as per customer requirement.
- Base and center column are made of Stainless Steel or Fibre Glass for easy cleaning and hygienic.
- Table frame, Side rails and Accessories are made up of stainless steel.

ACCESSORIES OF O.T. TABLE :

- Anaesthetic Frame (1Pcs)
- Padded Side Supports (1 Pair)
- Intravenous Armboard (1 Pair)
- Padded Shoulder Supports (1 Pair)
- Knee Crutches (1 Pair)
- One Matters Set Thickness 50mm



THE NEW BENCHMARK IN HEALTHCARE

VARIOUS POSITIONS OF TYCO 8004 A



Trendelenburg



Lateral Tilt (Both)



Reverse Trendelenburg



Flex Position



Chair Position



Reflex Position



Lithotomy Position



Kidney Position



Back Position

MINOR HEAD END CONTROL HYDRAULIC O.T. TABLE



TYCO 8002

TECHNICAL DATA : Minimum Height : 75 ± 7.5 cm
Maximum Height : 115 ± 5
Length Table Top : 180 cm to 195 cm
Width of Table Top : 50 cm
Lateral tilt : $20^\circ \pm 2^\circ$
Trendelenburg : $45^\circ \pm 3^\circ$
Chair Position : $-22^\circ / +90^\circ$
Flex/Reflex : $90^\circ / 22^\circ$
Kidney Position : 15 cm

SALIENT FEATURES : Single Action Hydraulic Lift
Heavy Mild Steel Base with Fibre wheels provides stability to the entire table
Auto floor locking device
The detachable leg section provides easy lithotomy Positioning
The improved Hydraulic system provides smooth movement and the speed of lowering the table is permanently set.

ACCESSORIES OF O.T. TABLE : Padded Side Supports (1 Pair)
Intravenous Armboard (1 Pair)
Padded Shoulder Supports (1 Pair)
Knee Crutches (1 Pair)
One Matters Set Thickness 50mm



THE NEW BENCHMARK IN HEALTHCARE

VARIOUS POSITIONS OF TYCO 8002



Trendelenburg



Chair Position



Reverse Trendelenburg



Lateral Tilt (Both)



Reflex



Flex

SIDE END CONTROL HYDRAULIC O.T. TABLE



TYCO 8001

TECHNICAL DATA : Minimum Height : 75 cm
Maximum Height : 95 cm
Length Table Top : 210 cm
Width of Table Top : 50 cm
Lateral tilt : $20^{\circ} \pm 2^{\circ}$
Trendelenburg : $30^{\circ} \pm 3^{\circ}$
Hydraulic Lift: 22 cm
Back Position: $+90^{\circ}/-40^{\circ}$
Kidney Position: 12 cm

SALIENT FEATURES : Single Action Hydraulic Lift.
Heavy Mild Steel Base with Fibre wheels.
Auto floor locking device.
The detachable leg section provides easy lithotomy Positioning.
The improved Hydraulic system provides smooth movement and the speed of lowering the table is permanently set.
Head and Leg section interchangeable.

ACCESSORIES OF O.T. TABLE : Anaesthetic Frame (1Pcs)
Padded Side Supports (1 Pair)
Intravenous Armboard (1 Pair)
Padded Shoulder Supports (1 Pair)
Knee Crutches (1 Pair)
One Matters Set Thickness 50mm



THE NEW BENCHMARK IN HEALTHCARE

VARIOUS POSITIONS OF TYCO 8001



Trendelenburg



Reverse Trendelenburg



Lateral Tilt (Both)



Chair Position



Flex / Reflex Position



Kidney Position

3600 REVOLVING HYDRAULIC O.T. TABLE as per ISI - 6106



TYCO 8001 R

TECHNICAL DATA :

- Minimum Height : 75 ± 5 cm
- Maximum Height : 95 ± 5 cm
- Length Table Top : 180 cm to 195 cm
- Width of Table Top : 50 cm
- Lateral tilt : $20^\circ \pm 2^\circ$
- Trendelenburg : $30^\circ \pm 3^\circ$
- Hydraulic Lift: 22 cm

SALIENT FEATURES :

- Single Action Hydraulic Lift.
- Heavy Mild Steel Base with Fibre wheels provides stability to the entire table .
- Auto floor locking device.
- The detachable leg section provides easy lithotomy Positioning.
- The improved Hydraulic system provides smooth movement and the speed of lowering the table is permanently set.

ACCESSORIES OF O.T. TABLE :

- Anaesthetic Frame (1Pcs)
- Padded Side Supports (1 Pair)
- Intravenous Armboard (1 Pair)
- Padded Shoulder Supports (1 Pair)
- Knee Crutches (1 Pair)
- One Matters Set Thickness 50mm



THE NEW BENCHMARK IN HEALTHCARE

VARIOUS POSITIONS OF TYCO 8001 R



Trendelenburg



Reverse Trendelenburg



Lateral Tilt (Both)



Chair Position

STANDARD ACCESSORIES



Shoulder Support



Anesthetic Screen Rod



Arm Board with Pad



Lateral Support



Knee Crutches Goepel type

ORTHOPEADIC ATTACHMENT OPTIONAL



OPTIONAL ATTACHMENTS



Body Elevator




Accessories Stand




Finger Traction



Knee Chest Position for Spine




Hand Traction Unit



Neck Elevator




Three Pin Fixture



Lateral Cassette Holder



U Shape Neuro Attachment



Shoulder Elevator